Treat each solution as follows. To $5.0 \, \text{ml}$ add $4 \, \text{ml}$ of dilute sulphuric acid R, $1 \, \text{ml}$ of ammonium molybdate solution R, $5 \, \text{ml}$ of water R and $2 \, \text{ml}$ of a solution containing $0.10 \, \text{g}$ of 4-methylaminophenol sulphate R, $0.5 \, \text{g}$ of anhydrous sodium sulphite R and $20.0 \, \text{g}$ of sodium metabisulphite R in $100 \, \text{ml}$ of water R. Shake and allow to stand for $15 \, \text{min}$. Dilute to $25.0 \, \text{ml}$ with water R and allow to stand for a further $15 \, \text{min}$. Measure the absorbance (2.2.25) at $730 \, \text{nm}$. Calculate the content of soluble phosphates from a calibration curve prepared using reference solutions (a), (b) and (c) after treatment.

Sulphates (2.4.13): maximum 0.6 per cent.

Dilute 8 ml of solution S to 100 ml with *distilled water R*. 15 ml of the solution complies with the limit test for sulphates.

Arsenic (2.4.2): maximum 1 ppm. 1.0 g complies with limit test A.

Heavy metals (2.4.8): maximum 20 ppm.

Dissolve 1.0 g in *dilute hydrochloric acid R* and dilute to 20 ml with the same acid. 12 ml of the solution complies with limit test A. Prepare the standard using *lead standard solution (1 ppm Pb) R*.

Loss on ignition. 10.0 per cent to 20.0 per cent, determined on 1.000 g at 800 ± 50 °C.

Neutralising capacity. Add 0.50 g to 30 ml of $0.1\,M$ hydrochloric acid previously heated to 37 °C and maintain at this temperature for 15 min while stirring. The pH (2.2.3) of the mixture after 15 min at 37 °C is 2.0 to 2.5.

ASSAY

Dissolve 0.400 g in 10 ml of dilute hydrochloric acid R and dilute to 100.0 ml with water R. To 10.0 ml of the solution, add 10.0 ml of 0.1 M sodium edetate and 30 ml of a mixture of equal volumes of ammonium acetate solution R and dilute acetic acid R. Boil for 3 min, then cool. Add 25 ml of alcohol R and 1 ml of a freshly prepared 0.25 g/l solution of dithizone R in alcohol R. Titrate the excess of sodium edetate with 0.1 M zinc sulphate until the colour changes to pink.

1 ml of 0.1 M sodium edetate is equivalent to 12.20 mg of AlPO $_{\scriptscriptstyle A}$.

STORAGE

In an airtight container.

IDENTIFICATION

- A. Solution S (see Tests) gives reaction (a) of sulphates (2.3.1).
- B. Solution S gives the reaction of aluminium (2.3.1).

TESTS

Solution S. Dissolve 2.5 g in *water R* and dilute to 50 ml with the same solvent.

Appearance of solution. Solution S is not more opalescent than reference suspension III (2.2.1) and is colourless (2.2.2, Method II).

pH (2.2.3): 2.5 to 4.0.

Dissolve 0.5~g in *carbon dioxide-free water R* and dilute to 25~ml with the same solvent.

Alkali and alkaline-earth metals: maximum 0.4 per cent.

To 20 ml of solution S add 100 ml of *water R*, heat and add 0.1 ml of *methyl red solution R*. Add *dilute ammonia R1* until the colour of the indicator changes to yellow. Dilute to 150 ml with *water R*, heat to boiling and filter. Evaporate 75 ml of the filtrate to dryness on a water-bath and ignite. The residue weighs a maximum of 2 mg.

Ammonium (2.4.1): maximum 500 ppm.

Dilute 0.4 ml of solution S to 14 ml with water R.

Iron (2.4.9): maximum 100 ppm.

Dilute 2 ml of solution S to 10 ml with *water R*. Use 0.3 ml of *thioglycollic acid R* in this test.

Heavy metals (2.4.8): maximum 50 ppm.

Dilute 8 ml of solution S to 20 ml with *water R*. 12 ml of the solution complies with test A. Prepare the reference solution using *lead standard solution* (1 ppm Pb) R.

ASSAY

Dissolve 0.500 g in 20 ml of *water R*. Carry out the complexometric titration of aluminium (2.5.11).

1 ml of 0.1 M sodium edetate is equivalent to 17.11 mg of $Al_2(SO_4)_3$.

STORAGE

In an airtight container.

01/2008:2156

01/2008:0165

ALUMINIUM SULPHATE

Aluminii sulfas

 $Al_2(SO_4)_3$, xH_2O M_r 342.1 (anhydrous substance)

DEFINITION

Content: 51.0 per cent to 59.0 per cent of $Al_2(SO_4)_3$. It contains a variable quantity of water of crystallisation.

CHARACTERS

 $\label{lem:appearance:colourless} Appearance : \ colourless, \ lustrous \ crystalls \ or \ crystalline \ masses.$

Solubility: soluble in cold water, freely soluble in hot water, practically insoluble in ethanol (96 per cent).

ALVERINE CITRATE

BY BRINE CITRIII

Alverini citras

$$CH_3$$
 CO_2H CO_2H

 $\substack{C_{26}H_{35}NO_7\\[5560-59-8]}$

DEFINITION

N-Ethyl-3-phenyl-*N*-(3-phenylpropyl)propan-1-amine dihydrogen 2-hydroxypropane-1,2,3-tricarboxylate.

Content: 99.0 per cent to 101.0 per cent (dried substance).

CHARACTERS

Appearance: white or almost white, crystalline powder. Solubility: slightly soluble in water and in methylene chloride, sparingly soluble in ethanol (96 per cent).

mp: about 104 °C.

IDENTIFICATION

Infrared absorption spectrophotometry (2.2.24).

Comparison: alverine citrate CRS.

TESTS

pH (2.2.3): 3.5 to 4.5.

Dissolve 0.250 g in *carbon dioxide-free water R* and dilute to 50.0 ml with the same solvent.

Related substances. Gas chromatography (2.2.28): use the normalisation procedure. *Use freshly prepared solutions*.

Test solution. Dissolve 0.250 g of the substance to be examined in water R and dilute to 20 ml with the same solvent. Add 2 ml of concentrated ammonia R and shake with 3 quantities, each of 15 ml, of methylene chloride R. To the combined lower layers add anhydrous sodium sulphate R, shake, filter, and evaporate the filtrate at a temperature not exceeding 30 °C, using a rotary evaporator. Take up the residue with methylene chloride R and dilute to 10.0 ml with the same solvent.

Reference solution (a). Dissolve 5 mg of alverine impurity D CRS (impurity D citrate) in 5 ml of water R, add 1 ml of concentrated ammonia R and shake with 3 quantities, each of 5 ml, of methylene chloride R. To the combined lower layers add anhydrous sodium sulphate R, shake, filter, and evaporate the filtrate at a temperature not exceeding 30 °C, using a rotary evaporator. Take up the residue with methylene chloride R, add 0.2 ml of the test solution and dilute to 2 ml with methylene chloride R.

Reference solution (b). Dilute 1.0 ml of the test solution to 100.0 ml with *methylene chloride R*. Dilute 1.0 ml of this solution to 20.0 ml with *methylene chloride R*.

Reference solution (c). Dissolve the contents of a vial of alverine for peak identification CRS (containing impurities C and E) in 1 ml of methylene chloride R.

Column:

material: fused silica;

- size: l = 25 m, $\emptyset = 0.32$ mm;

 stationary phase: poly(dimethyl)(diphenyl)siloxane R (film thickness 0.45 µm).

Carrier gas: helium for chromatography R.

Flow rate: 2.2 ml/min. Split ratio: 1:11. Temperature:

	Time	Temperature
Column	(min) 0 - 7	(°C) 120
Column	0 - 7	120
	7 - 13	$120 \rightarrow 240$
	13 - 21	240
	21 - 24	$240 \rightarrow 290$
	24 - 39	290
Injection port		290
Detector		290

Detection: flame ionisation.

Injection: 1 µl.

Identification of impurities: use the chromatogram supplied with *alverine for peak identification CRS* and the chromatogram obtained with reference solution (c) to identify the peaks due to impurities C and E.

Relative retention with reference to alverine (retention time = about 16 min): impurity A = about 0.28; impurity B = about 0.29; impurity C = about 0.46; impurity D = about 0.97; impurity E = about 1.7.

System suitability: reference solution (a):

 resolution: minimum 3.0 between the peaks due to impurity D and alverine.

Limits:

- impurities A, B: for each impurity, maximum 0.1 per cent;
- *impurity C*: maximum 0.2 per cent;
- *impurities D, E*: for each impurity, maximum 0.3 per cent;
- unspecified impurities: for each impurity, maximum 0.10 per cent;
- total: maximum 1.0 per cent;
- disregard limit: the area of the principal peak in the chromatogram obtained with reference solution (b) (0.05 per cent).

Heavy metals (2.4.8): maximum 20 ppm.

0.5 g complies with test G. Prepare the reference solution using 1 ml of *lead standard solution (10 ppm Pb) R*.

Loss on drying (2.2.32): maximum 0.5 per cent, determined on 1.000 g by drying in an oven at 80 $^{\circ}$ C for 2 h.

Sulphated ash (2.4.14): maximum 0.1 per cent, determined on 1.0 g.

ASSAY

Dissolve 0.375 g in 50 ml of *anhydrous acetic acid R*. Titrate with 0.1 M perchloric acid, determining the end-point potentiometrically (2.2.20).


1 ml of 0.1 M perchloric acid is equivalent to 47.36 mg of $C_{26}H_{35}NO_7$.

STORAGE

Protected from light.

IMPURITIES

Specified impurities: A, B, C, D, E.




A. R = Cl: 1-chloro-3-phenylpropane,

B. R = OH: 3-phenylpropan-1-ol,

C. $R = NH-C_2H_5$: N-ethyl-3-phenylpropan-1-amine,

D. N-(3-cyclohexylpropyl)-N-ethyl-3-phenylpropan-1-amine,



E. 3-phenyl-*N*,*N*-bis(3-phenylpropyl)propan-1-amine.